Instructions for Use

MRX DDU Calibrator, K5012

Contents

English [EN] _____	3
Čeština [CS] _____	6
Italiano [IT] _____	9
Slovenščina [SL] _____	12
Svenska [SV] _____	15

Instructions for Use [EN]

MRX D-dimer DDU Calibrator

REF K5012

For *In vitro* Diagnostic Use.

1 Intended use

For calibration of the MRX Red D-dimer, MRX Green D-dimer and MRX Blue D-dimer reagents in ng/mL D-dimer Units (DDU). Intended to be used by professional laboratory personnel using coagulation analysers.

2 Background and principle of method

For background information and principle of the method, refer to the Instructions for Use associated with each individual assay (K5011, K5034, K5035).

3 Components

MRX D-dimer DDU Calibrator consists of:
4 × 1 mL lyophilised citrated human plasma enriched with D-dimer.

4 Metrological traceability

Each new lot of MRX D-dimer DDU Calibrator is assigned with MRX Red D-dimer against in-house reference material with traceability to a working calibrator assigned according to ISO 17511:2003, section 5.6.¹ Refer to the Certificate of Analysis or the label for the lot-specific D-dimer concentration.

Specification:

Analyte	Specification
D-dimer DDU	2700–3100 ng/mL

5 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The calibrator contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

6 Preparation

- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water (e.g. MRX Laboratory Water, K5036). The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

7 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 24 h at 2 - 25 °C in the closed original vial, provided no contamination occurs.

8 Specimen collection and preparation

For specimen collection and preparation, refer to the Instructions for Use associated with each individual assay (K5011, K5034, K5035).

9 Procedure

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

10 Material required but not provided

Coagulation analyser capable of turbidimetric detection, pipettes and the following:

Reagent	REF
MRX Red D-dimer	K5034
MRX Green D-dimer	K5011
MRX Blue D-dimer	K5035

Control material	REF
MRX Routine Normal Control	K5039
MRX Routine Abnormal Control	K5040

Solutions	REF
Deionised water for reconstitution e.g. MRX Laboratory Water	K5036
Phosphate buffered saline (PBS) for dilution, e.g. MRX PBS Diluent	K5047

11 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Routine Controls (K5039/K5040) are recommended for MRX Red D-dimer, MRX Green D-dimer and MRX Blue D-dimer. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. Recalibration is suggested, as a minimum, whenever control plasmas are not within the acceptable range and each time a new batch of reagent is used.

12 Results

The results are reported in ng/mL D-dimer Units (DDU).

13 Expected values

For expected values, refer to the Instructions for Use associated with each individual assay (K5011, K5034, K5035).

14 Limitations and interfering substances

For limitations and interfering substances, refer to the Instructions for Use associated with each individual assay (K5011, K5034, K5035).

15 Analytical performance characteristics

For analytical performance characteristics, refer to the Instructions for Use associated with each individual assay (K5011, K5034, K5035).

16 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

17 Additional Information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

18 References

1. EN ISO 17511:2003 In vitro diagnostic medical devices - Measurement of quantities in biological samples - Metrological traceability of values assigned to calibrators and control materials

19 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains human blood or plasma derivatives

20 Revision history

Version	Changes to previous version
2.0	Addition of MRX Green D-dimer (K5011). Addition of languages.

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