

May-Grünwald solution

REF. 320070

Fixation and differential staining of cellular structures



IFU016B

Changes tracking 9

Legal representatives 9

For professional use only.

Read all information carefully before using this device.

IFU content may change, make sure you have the latest version available at

my.ral-diagnostics.fr.

Table of contents

Intended use.....	1
Principle.....	1
Device description	2
Storage and use conditions	2
Active components	2
Hazard classification and safety information.....	2
Personnel qualification	2
Specific equipment and reagents required but not provided	3
Operating procedure.....	3
Expected results.....	5
Performance	6
User quality control	6
Other products.....	6
Recommendations, notes, and troubleshooting.....	6
Table of symbols and abbreviations	8
Bibliography.....	8

Intended use

May-Grünwald solution is intended to be used in combination with Giemsa solution for the fixation and the differential staining of cellular structures prior microscopic examination.

If applicable, CellaVision RAL Diagnostics recommends using the associated CellaVision RAL Diagnostics products and cannot guarantee that the expected results will be achieved if used in combination with products of other brands.

Principle

Staining according to Pappenheim allows differential medullary and blood cells counting. It combines two stains: May-Grünwald and Giemsa stains. These are neutral mixtures with very distinct properties. They are not active in alcoholic medium and only act selectively when released in a buffered aqueous solution. This releasing induces the precipitation of neutral stains. May-Grünwald stains the acidophilic elements and the neutrophilic granulations of the leukocytes. Giemsa stains the cytoplasm of monocytes and lymphocytes as well as the chromatin of the nuclei.

Device description

May-Grünwald solution

Clear dark blue solution

REF. 320070-1000

1 x 1.0 L

REF. 320070-2500

1 x 2.5 L

For a specific batch, refer to the analysis certificate of the batch available at my.ral-diagnostics.fr.


Storage and use conditions

Storage and use temperature: 15-25°C.

Storage and use conditions: away from light and heat sources.

Bottle shelf life before opening: refer to the expiry date on the label.

Bottle shelf life after opening: refer to the expiry date on the label and if the "period after opening" symbol is present take it into account.



Active components

May-Grünwald solution

May- Grünwald: ca 0.3%

Hazard classification and safety information

May-Grünwald solution

Danger:

H225 - Highly flammable liquid and vapour.

H301+H311+H331 - Toxic if swallowed, in contact with skin or if inhaled.

H370 - Causes damage to organs.

P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P261 - Avoid breathing dust/fume/gas/mist/vapours/spray.

P264 - Wash hands thoroughly after handling.

P280 - Wear protective gloves, protective clothing, eye protection, face protection.

P301+P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor.

P308+P311 - IF exposed or concerned: Call a POISON CENTER or doctor.

CONT CH3OH



Personnel qualification

All samples and products must be handled by qualified and authorized personnel, using individual or collective protection, in accordance with the national directives in force in the laboratories. Personnel must also be aware of the classification of hazardous materials indicated on the label and the safety data sheet (available at my.ral-diagnostics.fr).

The diagnosis must be conducted by qualified and authorized personnel, in accordance with the procedures in force within the laboratory.

Specific equipment and reagents required but not provided

Sodium hyposulfite aqueous solution, ethanol, acetone, mounting media, microscope slides, and these following CellaVision RAL Diagnostics devices:

Giemsa R solution REF. 320310

pH=6.8 buffer solution for Hematology REF. 330368

pH=7.0 buffer solution for Hematology REF. 330370

Lugol, PVP-stabilized solution REF. 367400

This equipment may vary depending on the protocol. Refer to the relevant protocol (see the section operating procedure) to ensure that you have the necessary equipment to carry out tests.

Operating procedure

The equipment used for sample processing must comply with the Instructions for use of the supplier.

Sample preparation

The specimen must be treated in accordance with procedures available in the laboratory and promulgated by national authorities.

Manual blood smear: Mix the tube by slow inversion and install a smearing droplet device. Invert the tube and lightly press the drop depositor onto a slide to deposit a small drop of blood (Fig. 1-slide A at step 1).

Using another slide tilted at 45° (Fig. 1- slide B at step 1), spread the blood by capillarity on the short edge (Fig. 1-steps 2 and 3) using a pushing motion (Fig. 1-step 4). A good quality smear does not reach the end of the slide and has a gradual decrease in thickness until the end is feathered. Allow the smear to air dry before fixing or staining.

Note: if you do not have a smearing droplet device, open the tube, and use a pipette to deposit a blood drop.




Figure 1. Schematic representation of performing a blood smear

A & B: slides, 1-4: steps 1 to 4

Manual bone marrow smear by crushing method: using a pipette, deposit a small amount of the sample on a microscope slide. Blot up excess blood to keep only shiny lumps. Cover the first slide with a slide. Squeeze and thin the sample by sliding and stretching to the end of the slide. A good quality smear does not reach the end of the slide. Discard the slide used for smearing. Allow the smear to air dry before fixing or staining.

Histological sections: dewax and et hydrate tissues sections in appropriate reagents before staining.

Reagents and instruments preparation


If applicable dilute the May-Grünwald solution and the Giemsa solution according to the indications in the protocol section. Transfer the solutions into staining baths as indicated in the protocols below.

Acetic water: 5 drops in 100 mL of distilled water

5% sodium hyposulfite aqueous solution: 5 g of sodium hyposulfite in 100 mL of distilled water

Protocols

The staining steps of the protocols indicated below consist of a successive covering of the slides with the different staining reagents or dipping of the slides in the different staining baths. Refer to the title to know which case you are in. The processing time only considers the contact time with the reagents.



Note:

For bath staining method, fill each jar with the appropriate reagent in enough volume to cover the smears. Be sure to cover the entire sample when immersing the slide.



Note:

For the covering staining method, place the slide on a stand with the fixed smear on top.

Protocol hematological samples - Manual bath staining method - Manual microscopic analysis

Processing time [hh: mm: ss]: 00: 14 :10

Steps	Reagent	Time [mm:ss]	Indications
Fix and prestain	May-Grünwald solution	03:00	
Rinse	Buffer Solution	01:00	N/A
Stain	1/20 Giemsa R solution diluted in buffer solution	10:00	
Rinse	Buffer Solution	00:10	Agitate continuously in the bath during the countdown.
Dry	N/A	≥03:00	N/A

Protocol hematological samples - Manual covering staining method - Manual microscopic analysis

Processing time [hh: mm: ss]: 00: 14 :10

Steps	Reagent	Time [mm:ss]	Indications
Fix and prestain	1 mL of May-Grünwald solution	03:00	Cover the slide
Rinse	1 mL of Buffer Solution	01:00	Add carefully and mix without spilling. Get rid of the reagent at the end of countdown.
Stain	1/30 Giemsa R solution diluted in buffer solution	10:00	Cover the slide and get rid of the reagent at the end of the countdown.
Rinse	Buffer Solution	00:10	Quickly. Tap water can be used.
Dry	N/A	≥03:00	N/A

Protocol for histological sections - Manual bath staining method- Manual microscopic analysis

Processing time [hh: mm: ss]: 01: 05 :00

Steps	Reagent	Time [mm: ss]	Indications
Stain	Lugol, PVP-stabilized solution	05:00	N/A
Stain	5% sodium hyposulfite aqueous solution	05:00	N/A
Rinse	Distilled water	N/A	
Stain	1/5 May-Grünwald solution diluted in distilled water	15:00	In an autoclave at 37° C
Stain	Add Giemsa L solution (3 drops in 2 mL of distilled water)	40:00	
Rinse	Distilled water	N/A	N/A
Differentiate	Acetic water	N/A	N/A
Rinse	Distilled water	N/A	Then drain the excess onto filter paper
Dehydrate	Ethanol/ acetone mixture 50/50	N/A	N/A
Dehydrate	Toluene or xylene	N/A	2 baths
Mount	Toluene or xylene-based mounting media	N/A	N/A

Expected results

Hematological samples

Nuclei/chromatin: +/- dense purple
Granulocytes cytoplasm without RNA: light purplish-pink
Granulocytes eosinophilic granules: orangey-pink
Granulocytes basophilic granules: dark blue
Granulocytes neutrophilic granules: ± deep purple-pink
Lymphocytes cytoplasm with RNA: pure blue
Lymphocytes cytoplasm without RNA: light blue
Lymphocytes azurophilic granules: red
Monocytes cytoplasm: purplish-blue
Erythrocytes: pinkish-beige to beige-grey
Platelets chromomere: purplish-red
Platelets hyalomere: bluish
Blood parasites nucleus (*Plasmodium*): red
Blood parasites cytoplasm (*Plasmodium*): blue

Histological sections samples

Nuclei/Chromatin: purple to pink
Basophilic cytoplasm: sky to dark blue
Acidophilic cytoplasm: light red to pinky
Polychromatophilic cytoplasm: greyish or purplish
Acidophilic leukocyte granulations: orangey
Neutrophilic leukocyte granulations: dirty brown pink
Basophilic leukocyte granulations: dark violet
Azurophilic leukocyte granulations: purple or purplish
Basophilic erythrocyte granulations: cobalt blue

If observed results vary from those expected, contact CellaVision RAL Diagnostics technical service through your usual supplier for assistance.

Performance

May-Grünwald solution enables cell structure staining and microscopic analysis. As it does not allow the detection of analytes, analytical performance is not applicable to this reagent.

This medical device is based on scientific validity (scientific peer-reviewed literature) and demonstration of clinical performance through experience gained from routine diagnostic testing, as well as the regular evaluation of these performances within the framework of Post Market Performance Follow-up (PMPF), to ensure that it continues to meet the expected performance and safety standards.

To ensure product performance, use clean and dry laboratory equipment.

The laboratory is responsible for notifying the manufacturer and state competent authority of any serious incident relating to the medical device uses.

User quality control

Users are responsible for determining the appropriate quality control procedures for their laboratory and complying with applicable laboratory regulations.

Hematological sample: CellaVision RAL Diagnostics recommends staining a freshly made blood smear with a normal WBC count and no known abnormal pathology at reagent renewal and for the first staining cycle each day. Slides stained for quality control purposes should be checked to ensure that they are satisfactory for the intended test (properly stained and free of precipitates).

Staining results must also be compliant with this manual expected results.

These quality control procedures should only be performed by qualified personnel.

Other products

For more information, contact your usual supplier.

Recommendations, notes, and troubleshooting

Products appearance

If the appearance of the products differs from the description above, do not use it and contact CellaVision RAL Diagnostics technical service through your usual supplier for assistance.

Procedures notes

To prevent products degradation, comply with the storage and handling recommendations specified in this manual.

Be sure to cover the entire sample when immersing the slide.

Staining quality and reproducibility are obtained by using a buffer solution.

Due to the complex mode of action of stains, it is essential to set up standard staining conditions to ensure perfect staining quality and reproducibility. The use of tap water or water mixing is not recommended because unpredictable variations may occur and alter staining results.

To ensure staining quality and reproducibility, CellaVision RAL Diagnostics recommends the use of a buffer solution specially formulated for Hematology. Safe for the users, Buffer solutions for Hematology are formulated with phosphates, are specially developed for hematology, and guarantee improved product stability and rinsability.

Products stability

Every CellaVision RAL Diagnostics product can be used until the expiry date indicated on, in its original packaging if it is still hermetically sealed.

Staining stability

Staining quality and reproducibility depend on the correct use of the products. Staining conducted according to these recommendations remain stable for several days. If it is necessary to store the stained smears for several months or years, CellaVision RAL Diagnostics recommended mounting them with a coverslip, using a suitable mounting liquid and storing them in a light and dustproof container.

Instructions for cleaning and waste disposal

All biological samples, effluents, and used consumables should be considered potentially hazardous.





































To avoid any risk, apply the following instructions: dispose of samples, effluents, and consumables in accordance with laboratory standards and applicable national and local standards and regulations.

Chemical and biological waste must be collected and processed by specialized, registered companies.

Table of symbols and abbreviations

Depending on the product, you may find the following symbols on the device or the packaging material.

GHS Pictograms	Interpretation
	Explosive
	Flammable
	Oxidizer
	Compressed gas
	Corrosive
	Toxic
	Harmful
	Health Hazard
	Environmental Hazard
	No labelling applicable

SYMBOL	INTERPRETATION
	Batch code
	Serial number
	Catalogue reference
	Date of manufacture
	Use up to
	Unique device identifier
	Manufacturer
	Importer
	Entity distributing the medical advice in the region concerned
	CE marking device
	In vitro diagnostic medical device
	Authorised Representative in the European Community
	Authorised Representative in Switzerland
	Authorised Representative in United Kingdom
	Complies with UK guidelines
	Do not use if packaging is damaged
	Keep away from light Keep away from heat
	Temperature limit: 15-25°
	Temperature limit: 15-30°
	Keep dry
	Box: handling upwards
	Fragile
	Sterilised by irradiation
	Single sterile barrier system with outer protective packaging
	Sterile and radiation-sterilised barrier suit
	Do not reuse
	Do not resterilize
	Contents sufficient for n tests
	Hazardous material contained
	Consult instructions for use
	Use
	After opening, use within XX months
	The product must not be used in conjunction with an automatic colouring machine
	Indicates a medical device that contains potentially carcinogenic, mutagenic or reprotoxic (CMR) substances, or substances classified as endocrine disruptors

Bibliography

1. **DUHAMEL G., DUHAMEL E.**, *Cytologie hématologique, Les cellules pathologiques I et II, Coloration au May-Grünwald Giemsa RAL, Biologiste et Praticien et Réactifs RAL*, 1984 et 1989.
2. **Ecole Nationale de Chimie**, *Coloration de Pappenheim, Présentation théorique des mécanismes cytochimiques des colorants neutres avec applications techniques détaillées, Journée du technicien biologiste*, mars 1980, p. 1-9.
3. **GANTER P., JOLLES G.**, *Histochimie normale et pathologique et pathologique*, éd GAUTHIER-VILLARS, vol 2, 1970 p. 1435-1436
4. **GENTILHOMME O., TREILLE-RITOUET D., BRYON P-A.**, *Cytologie hématologique, Les cellules normales, Coloration au May-Grünwald Giemsa RAL, Réactifs R.A.L.*, 1989.
5. **LANGERON M.**, *Précis de microscopie*, Masson & Cie 6^{ème} éd., 1942, 587-591.
6. **MATHIOT C.**, *Cytologie en hématologie, quelques aspects de la pathologie*. Biologiste, praticien et Réactifs RAL, 1979
7. **SOCIETE FRANCAISE D'HEMATOLOGIE (SFH)**, *Guides des bonnes pratiques des ponctions médullaires*, Juin 2003, VI.2
8. **THEML H.**, *ATLAS de poche d'Hématologie, Médecine-Sciences Flammarion*, p. 19-25, 2000

Changes tracking

Date	Version	Changes
12/2024	IFU016B	Update in header and the following paragraphs: Device description, Storage and use conditions, Hazard classification and safety information, Operating procedure, Expected results, Performance and Table of symbols and abbreviations. Add of the covering staining method. Add of the legal representatives. Add of CH-REP and UK-REP symbols.
05/2022	IFU016A	IVDR (EU) 2017/746 compliance

Legal representatives

Countries	Address
UK REP	Qavis UK Ltd, company N° SC679796, 56-66 Frederick Street Edinburgh, EH21LS, United Kingdom
CH REP	MedEnvoy Switzerland, Gotthardstrasse 28, 6302 Zug Switzerland

